Felivet Feline Immunodeficiency Virus Antibody – Feline Leukemia Virus Antigen Combo Test

■ INTENDED USE

Feline Immunodeficiency Virus Antibody – Feline Leukemia Virus Antigen Combo Test is a lateral flow immunochromatographic assay for the feline immunodeficiency virus antibody (FIV Ab) and for the qualitative detection of feline leukemia virus antigen (FeLV Ag) in cat's blood specimen.

Assay Time: 5-10 minutes

■ PRINCIPLE

The Feline Immunodeficiency Virus Antibody – Feline Leukemia Virus Antigen Combo Test is based on sandwich method lateral flow immunochromatographic assay.

■ REAGENTS AND MATERIALS

- 5 test pouches, with cards and disposable droppers
- 1 bottle of assay buffer
- 5 capillary droppers
- 1 package insert

■ STORAGE AND STABILITY

The kit can be stored at room temperature (2-30°C). The test kit is stable through the expiration date (24 months) marked on the package label. **DO NOT FREEZE.** Do not store the test kit in direct sunlight.

■ SPECIMEN PREPARATION AND STORAGE

1. Specimen should be obtained and treated as below.


Serum or plasma: collect the whole blood for the patient cat, centrifuge it to get the serum, or place the whole blood into a tube which contains anticoagulants to get plasma.

Whole blood: directly use the whole blood in the assay.

2. All specimen should be tested immediately. If not for testing right now, they should be stored at 2-8°C.

■ TEST PROCEDURE

- Allow all materials, including specimen and test device, recover to 15-25°C before running the assay.
- Take out the test card from the foil pouch and place it horizontally.
- Using the capillary dropper to collect the prepared specimen (serum, plasma or whole blood), and place 1 drop (10 µL) of the into the sample well "S" of the test card. Then place 3 drops (approx. 90 µL) of assay into the sample well. Start the timer.
- Interpret the result in 5-10 minutes. Result after 10 minutes is considered as invalid.

■ INTERPRETATION OF RESULTS

- Positive (+): The presence of both "C" line and zone "T" line, no matter T line is clear or vague.
- Negative (-): Only clear C line appear. No T line.
- Invalid: No colored line appears in C zone. No matter if T line appears.

■ PRECAUTIONS

- All reagents must be at room temperature before running the assay.
- Do not remove test cassette from its pouch until immediately before use.
- Do not use the test beyond its expiration date.
- The components in this kit have been quality control tested as standard batch unit. Do not mix components from different lot numbers.
- All specimens are of potential infection. It must be strictly treated according to the rules and regulations by local states.

■ LIMITATION

Feline Immunodeficiency Virus Antibody – Feline Leukemia Virus Antigen Combo Test is designed for in vitro veterinary diagnosis use only. All result should be considered with other clinical information available with veterinarian. It is suggested to apply a further confirmative method such as Western Blot when positive result was observed.

Felivet Feline Immunodeficiency Virus Antibody – Feline Leukemia Virus Antigen Combo Test

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